



...in an emergency

Evacuation Aids

- First Aid
- Fire, Health & Safety Aids



ESCAPE

MOBILITY COMPANY

a trademark of van Leeuwen Companies

Horizontal evacuation

First Aid

Our range of 'quality proof' first aid kits have been chosen to give you peace of mind so that you can be assured that your business is fully compliant with all of the latest Health & Safety Regulations.

Evacuation Aids

Escape Mobility Company are the specialists in evacuation aids amongst which includes the Escape-Chair®, an indispensable solution for transporting people down stairs in an emergency. Now added to our range is the introduction of the Escape-Mattress and Escape Sheet.

Fire, Health & Safety Aids

Our range of 'quality proof' fire, health & safety aids have been chosen as being the most sought after products that our customers have identified, following completion of their Fire Risk Assessments. So we bring you a small range of these products to make your life easier and to help you comply to the Regulatory Reform (Fire Safety) Order 2005.

P.O. box, NL-6401 DD Heerlen
T + 44 (0)207 2450695
F + 44 (0)207 2450085
E info@escape-mobility.com
W www.escape-mobility.com

ESCAPE

MOBILITY COMPANY

a trademark of van Leeuwen Companies



QUALITY
PROOF
MOBILITY COMPANY

ARC2007/6341UK



Horizontal evacuation ...in an emergency

Evacuation Aids

Being a specialist in evacuation aids and our branded Escape-Chairs, we have extended our range of products to include horizontal transport aids for bringing mobility-impaired people down the stairs to a place of safety. For this purpose we have developed the Escape-Mattress and Escape-Sheet as well as their corresponding storage facilities.

Quality Proof of the Mobility Company

As well as our own products, we offer a range of other products that have been tested thoroughly to ensure they meet our high quality standards. That is why they have the 'Quality Proof' mark of the Mobility Company.



ESCAPE

MOBILITY COMPANY

a trademark of van Leeuwen Companies

ESCAPE-MATTRESS FOLD

Technical Details:

- Measurements (folded)
 - Height: ca. 97 cm
 - Width: ca. 68 cm
 - Thickness: ca. 15 cm
- Measurements (unfolded)
 - Length: ca. 193 cm
 - Width: ca. 68 cm
 - Thickness: ca. 7 cm
- Detachable and adjustable safety head rest
- Incl. webbed strapping

ESCAPE-LOCKER

- Secure storage especially designed for your evacuation chair, mattresses or sheets. Has a secure lock and is supplied with two keys.
- Can be fixed to the wall.
- Obtainable in white and green
- Measurements: 130 x 80 x 30 cm



ESCAPE-MATTRESS COMPACT

Technical Details:

- Measurements (folded)
 - Height: ca. 65 cm
 - Width: ca. 65 cm
 - Thickness: ca. 15 cm
- Measurements (unfolded)
 - Length: ca. 200 cm
 - Width: ca. 65 cm
 - Thickness: ca. 7 cm
- Detachable and adjustable safety head rest
- Incl. webbed strapping
- Measurements storage bag: ca. 73 x 83 x 25 cm



ESCAPE-SHEET

The Escape-Sheet is easily placed between the bed and the mattress and secured with elastic straps over the corners of the mattress. In an emergency, strong webbed safety belts are strapped over the top of the patient and bed-clothing, securing them to the mattress. The Escape-Sheet, together with mattress and patient, is then pulled easily onto the floor and down-stairs, to a place of safety. The Escape-Sheet is suitable for use in Hospitals, Care and Nursing Homes.



Technical details:

- Dimensions: 220 x 90 cm
- Made from 100% 3-mesh polyester
- Washable.



CARRY SHEET SPECIAL

Safety sheet (with footsack). Fitted with 8 hand grips. Ca. 200 x 70 cm



CARRY SHEET SPECIAL XXL

Safety sheet (with footsack). Fitted with 14 hand grips. Ca. 220 x 130 cm



STRETCHER

Foldable aluminium stretcher. Supplied with head rest. Ca. 104 x 56 x 14 cm (folded) Ca. 230 x 56 x 14 cm (unfolded)



MINI-STRETCHER

Flexible mini-stretcher that helps transfer people, who may be lying on the ground or in a bed.



We reserve the right to modify the technical and optical specifications.

